

# Tabla de revelado

## Argenti NANOTOMIC-X

Revelador	Dilución	Temp.	Tiempo en min.	ISO
Rodinal	1+200	20°	60	25
R09 Studio (Studional)	1+45	20°	7	12
Argenti Nanodol	1+19	25°	8 - 10	50
Argenti Nanodol	1+19	25°	10 - 12	100
Argenti Nanodol	1+19	25°	15 - 18	200

[www.foto-r3.com](http://www.foto-r3.com)



FOTO R3